

## Product datasheet for **MP202418**

### Cxcl2 Mouse qPCR Primer Pair (NM\_009140)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20310
Forward Sequence:	CATCCAGAGCTTGAGTGTGACG
Reverse Sequence:	GGCTTCAGGGTCAAGGCAAAC
ACCN:	<u><a href="#">BC119511</a></u> , <u><a href="#">BC119512</a></u> , <u><a href="#">NM_009140</a></u> , <u><a href="#">NM_009140.1</a></u> , <u><a href="#">NM_009140.2</a></u>
UniProt ID:	<u><a href="#">P10889</a></u>
Synonyms:	CINC-2a; Gro2; GROb; Mgsa-b; MIP-2; MIP-2a; Mip2; Scyb; Scyb2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)