

Product datasheet for **MP202425**

Ackr3 Mouse qPCR Primer Pair (NM_007722)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12778
Forward Sequence:	GACCGCTATCTCTCCATCACCT
Reverse Sequence:	GTTGGAAGCAGATGTGACCGTC
ACCN:	<u>BC015254</u> , <u>NM_007722</u> , <u>NM_007722.1</u> , <u>NM_007722.2</u> , <u>NM_007722.3</u> , <u>NM_007722.4</u> , <u>BQ033958</u>
UniProt ID:	<u>P56485</u>
Synonyms:	AW541270; Cmkor1; CXC-R7; CXCR-7; Cxcr7; RDC-1; Rdc1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)