

Product datasheet for **MP202427**

Cklf Mouse qPCR Primer Pair (NM_029295)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75458
Forward Sequence:	GCTCATTGTGTCTGTGTTGGCTC
Reverse Sequence:	CCGAAGTTTCTGGCACATAAGGG
ACCN:	BC027684 , NM_029295 , NM_029295.1 , NM_029295.2 , BE630757
UniProt ID:	Q9DAS1
Synonyms:	1700001C14Rik; 1810018M11Rik; A730028G07Rik; C32; CKLF1; CKlf2; CKLF3; CKLF4; CKLF5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)