

Product datasheet for **MP202444**

Clcn4 Mouse qPCR Primer Pair (NM_011334)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12727
Forward Sequence:	GGAGGTGTGCTCTTCAGTCTAG
Reverse Sequence:	AGAGAACCAGGCGCTATTTCC
ACCN:	<u>BC005553</u> , <u>BC110668</u> , <u>NM_011334</u> , <u>NM_011334.1</u> , <u>NM_011334.2</u> , <u>NM_011334.3</u> , <u>NM_011334.4</u> , <u>BC089497</u> , <u>BQ176890</u> , <u>BU610450</u> , <u>BY035169</u> , <u>BY790669</u>
UniProt ID:	<u>Q61418</u>
Synonyms:	cIC-4; Clc4-2; Clcn4-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)