

Product datasheet for **MP202466**

Cckbr Mouse qPCR Primer Pair (NM_007627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12426
Forward Sequence:	TCACACTCCTGCCCAACCTCAT
Reverse Sequence:	CCGCGAGATTTAGAGTGGACAC
ACCN:	BC103530 , BC103531 , BC103538 , NM_007627 , NM_007627.1 , NM_007627.2 , NM_007627.3 , NM_007627.4 , NM_007627.5 , BC103537 , BY245616
UniProt ID:	P56481
Synonyms:	CCK-BR; CCK2; CCK2-R; CCK2R; CCKR-; CCKR-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)