

Product datasheet for **MP202508**

Chrac1 Mouse qPCR Primer Pair (NM_053068)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93696
Forward Sequence:	AGGTGTCGAGCATCAACCAGGA
Reverse Sequence:	CTTGGCTTTGCCACTTCCGTGT
ACCN:	<u>BC061043</u> , <u>NM_053068</u> , <u>NM_053068.1</u> , <u>NM_053068.2</u> , <u>NM_053068.3</u> , <u>BC016593</u>
UniProt ID:	<u>Q9JKP8</u>
Synonyms:	2410152E03Rik; 2810406L04Rik; YCL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)