

## Product datasheet for **MP202511**

### Cdt1 Mouse qPCR Primer Pair (NM\_026014)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67177
Forward Sequence:	CTCCACAATCGCTCTGAGACTG
Reverse Sequence:	GGAAGCGATACGACGTGGGATA
ACCN:	<u><a href="#">BC048076</a>, <a href="#">NM_026014</a>, <a href="#">NM_026014.1</a>, <a href="#">NM_026014.2</a>, <a href="#">NM_026014.3</a>, <a href="#">BC013842</a>, <a href="#">BC024485</a></u>
UniProt ID:	<u><a href="#">Q8R4E9</a></u>
Synonyms:	2610318F11Rik; AW545653; C76791; Ris2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)