

Product datasheet for **MP202549**

Crocc Mouse qPCR Primer Pair (NM_172122)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	230872
Forward Sequence:	ATCGCCAAGAGGCTAGTGAGCT
Reverse Sequence:	TCTGCTCCTTGCTCCTCATTGGC
ACCN:	BC054054 , NM_172122 , NM_172122.1 , NM_172122.2 , BC027090 , BC034536 , BC036729 , BC038074 , BQ937768
UniProt ID:	Q8CJ40
Synonyms:	KIAA0445
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)