

Product datasheet for **MP202551**

Clock Mouse qPCR Primer Pair (NM_007715)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12753
Forward Sequence:	GGCTGAAAGACGGCGAGAACTT
Reverse Sequence:	GTGCTTCCTTGAGACTCACTGTG
ACCN:	<u>NM_007715</u> , <u>NM_007715.1</u> , <u>NM_007715.2</u> , <u>NM_007715.3</u> , <u>NM_007715.4</u> , <u>NM_007715.5</u> , <u>NM_007715.6</u> , <u>BC129886</u> , <u>BC148468</u> , <u>BC156741</u> , <u>BX526669</u>
UniProt ID:	<u>O08785</u>
Synonyms:	5330400M04Rik; KAT13D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)