

Product datasheet for **MP202571**

Clint1 Mouse qPCR Primer Pair (NM_001045520)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216705
Forward Sequence:	ACCTGGTATGCAGCCTTCCAAG
Reverse Sequence:	GGAGAAGTGGGTTACAGCTC
ACCN:	<u>NM_001045520</u> , <u>NM_001045520.1</u> , <u>NM_001045520.2</u> , <u>NM_001045520.3</u> , <u>BC004080</u> , <u>BC167227</u> , <u>BX514666</u> , <u>BY628486</u> , <u>NM_001045520.4</u>
Synonyms:	AI642036; AW537616; C530049I24Rik; Enth; Epn4; F630206G17Rik; mKIAA0171
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)