

Product datasheet for **MP202579**

Cldn13 Mouse qPCR Primer Pair (NM_020504)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57255
Forward Sequence:	AGATGCCACCATTGGGAGGGT
Reverse Sequence:	AGCAATCCCAGCCAGGTTCCAA
ACCN:	<u>BC115481</u> , <u>BC115482</u> , <u>NM_020504</u> , <u>NM_020504.1</u> , <u>NM_020504.2</u> , <u>NM_020504.3</u> , <u>NM_020504.4</u> , <u>BY079003</u>
UniProt ID:	<u>Q9Z0S4</u>
Synonyms:	2500004B14Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)