

Product datasheet for **MP202627**

F2r Mouse qPCR Primer Pair (NM_010169)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14062
Forward Sequence:	TTTCCTCTCCGACAGTCCTGGT
Reverse Sequence:	CAGCACAAAGATGCTGTAGAGGTG
ACCN:	<u>BC031516, NM_010169, NM_010169.1, NM_010169.2, NM_010169.3, NM_010169.4</u>
UniProt ID:	<u>P30558</u>
Synonyms:	AI482343; Cf2r; Par1; ThrR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)