

Product datasheet for **MP202651**

Coasy Mouse qPCR Primer Pair (NM_027896)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71743
Forward Sequence:	GGAAACCTGCTTCAGCTCCAA
Reverse Sequence:	ATGCTCCCAAGTTCTTCAGCCG
ACCN:	<u>NM_027896</u> , <u>NM_027896.1</u> , <u>NM_027896.2</u> , <u>NM_027896.3</u> , <u>NM_027896.4</u> , <u>BB872569</u> , <u>BC020046</u> , <u>BC064753</u> , <u>BC156119</u> , <u>BC157111</u>
Synonyms:	1300003G02Rik; Dpck; Ppat; Ukr1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)