

Product datasheet for **MP202671**

Ccdc102a Mouse qPCR Primer Pair (NM_001033533)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	234582
Forward Sequence:	CCAGAAAGTGCTCCTCAAGGAG
Reverse Sequence:	GAGCATCTCCTGCTTGGTCTTG
ACCN:	NM_001033533 , NM_001033533.1 , NM_001033533.2 , NM_001033533.3 , BC139307 , BB613877 , BC027663 , BC107374 , BC107375 , BC114992 , BC139308 , BF020220
UniProt ID:	Q3TMW1
Synonyms:	BC027663; Gm586
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)