

Product datasheet for **MP202719**

Ccdc34 Mouse qPCR Primer Pair (NM_026613)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68201
Forward Sequence:	AGGAATAGAGTTCCAGCAGCGC
Reverse Sequence:	CCTTGCCAACAAACCACATCTCC
ACCN:	<u>BC119347</u> , <u>BC120580</u> , <u>NM_026613</u> , <u>NM_026613.1</u> , <u>NM_026613.2</u> , <u>NM_026613.3</u> , <u>NM_026613.4</u> , <u>BC049229</u> , <u>BC062248</u> , <u>BC089458</u> , <u>BC108339</u> , <u>BQ748443</u>
UniProt ID:	<u>Q3UI66</u>
Synonyms:	2810027O19Rik; 4930522P10Rik; AI596202
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)