

Product datasheet for **MP202727**

Cep95 Mouse qPCR Primer Pair (NM_177088)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	320162
Forward Sequence:	CCAAACGGAGTCAGGTCCAGAA
Reverse Sequence:	GCATCAGGGATAGCAGACTGTG
ACCN:	<u>BC116308</u> , <u>NM_177088</u> , <u>NM_177088.1</u> , <u>NM_177088.2</u> , <u>NM_177088.3</u> , <u>BC090675</u> , <u>BC116309</u> , <u>BE852789</u> , <u>BU695983</u>
UniProt ID:	<u>Q8BVV7</u>
Synonyms:	4732496G21Rik; AI118346; AI448335; Ccdc45; D330027A14Rik; F630025I20Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)