

Product datasheet for **MP202793**

Csde1 Mouse qPCR Primer Pair (NM_144901)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	229663
Forward Sequence:	AAGAAGCAGGCTCGCTGTCAAG
Reverse Sequence:	GTCATCTCCAGGCTGTAGGGTT
ACCN:	BC016898 , BC029058 , BC029688 , BC031404 , BC062097 , NM_144901 , NM_144901.1 , NM_144901.2 , NM_144901.3 , NM_144901.4 , BC024826 , BC036720 , BC055072
UniProt ID:	Q91W50
Synonyms:	AA960392; BC016898; D3Jfr1; mKIAA0885; unr
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)