

Product datasheet for **MP202801**

Col3a1 Mouse qPCR Primer Pair (NM_009930)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12825
Forward Sequence:	GACCAAAAGGTGATGCTGGACAG
Reverse Sequence:	CAAGACCTCGTGCTCCAGTTAG
ACCN:	BC043089 , BC052398 , BC058724 , NM_009930 , NM_009930.1 , NM_009930.2 , BC013626 , BC028248 , BC039923 , BY139268
UniProt ID:	P08121
Synonyms:	AW550625; Col3; Col3a-1; Tsk; Tsk-; Tsk-2; Tsk2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)