

Product datasheet for **MP202842**

Csf1r Mouse qPCR Primer Pair (NM_001037859)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12978
Forward Sequence:	TGGATGCCTGTGAATGGCTCTG
Reverse Sequence:	GTGGGTGCATTCCAACCTGC
ACCN:	<u>BC036343</u> , <u>BC043054</u> , <u>BC066863</u> , <u>NM_001037859</u> , <u>NM_001037859.1</u> , <u>NM_001037859.2</u> , <u>BB738990</u> , <u>BC050024</u>
UniProt ID:	<u>P09581</u>
Synonyms:	AI323359; CD115; CSF-1R; Csfmr; Fim-2; Fim2; Fms; M-CSF-R; M-CSFR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)