

Product datasheet for **MP202857**

Commd8 Mouse qPCR Primer Pair (NM_178599)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27784
Forward Sequence:	GCTGAGTGAAGCATGTTCTCG
Reverse Sequence:	CAAAACAGCTTCTTGGTGACATGA
ACCN:	BC098183 , NM_178599 , NM_178599.1 , NM_178599.2 , NM_178599.3 , NM_178599.4 , BC023888 , BC051517 , NM_178599.5
UniProt ID:	Q9CZG3
Synonyms:	C77770; D5Buc26e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)