

## Product datasheet for **MP202871**

### C3 Mouse qPCR Primer Pair (NM\_009778)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 12266  |
| Forward Sequence: | CGCAACGAACAGGTGGAGATCA   |
| Reverse Sequence: | CTGGAAGTAGCGATTCTTGGCG   |
| ACCN:             | <u><a href="#">BC043338</a></u> , <u><a href="#">NM_009778</a></u> , <u><a href="#">NM_009778.1</a></u> , <u><a href="#">NM_009778.2</a></u> , <u><a href="#">NM_009778.3</a></u> , <u><a href="#">BC029976</a></u> , <u><a href="#">BY326535</a></u>  |
| UniProt ID:       | <u><a href="#">P01027</a></u>  |
| Synonyms:         | AI255234; ASP; HSE-MSF; P1p  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)