

Product datasheet for **MP202934**

Cpne8 Mouse qPCR Primer Pair (NM_025815)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66871
Forward Sequence:	ACGCAGATCAACTTCACGGTGG
Reverse Sequence:	TCTCCTACAGCCTTCAGTGCCA
ACCN:	BC017626 , NM_025815 , NM_025815.1 , BC017626.1 , BB043122 , BC048551 , BC076564 , BC138745 , BC138746 , BU846828
UniProt ID:	Q9DC53
Synonyms:	1200003E11Rik; 1500031E20Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)