

Product datasheet for **MP202962**

Crhr1 Mouse qPCR Primer Pair (NM_007762)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12921
Forward Sequence:	CGCAAGTGGATGTTCTGTGCA
Reverse Sequence:	TCCAGGACGTTTGCCAAACCAG
ACCN:	BC103675 , BC104735 , BC105673 , BC105675 , NM_007762 , NM_007762.1 , NM_007762.2 , NM_007762.3 , NM_007762.5 , BQ174858 , BY278495
UniProt ID:	P35347
Synonyms:	CRF; CRF-R; CRF1R; CRFR1; Crhr
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)