

## Product datasheet for **MP202979**

### Crebbp Mouse qPCR Primer Pair (NM\_001025432)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12914
Forward Sequence:	CACCATCTGTGGCTACTCCTCA
Reverse Sequence:	GGTTTCAGCACTGGTCACAGAG
ACCN:	<a href="#">NM_001025432</a> , <a href="#">NM_001025432.1</a> , <a href="#">BC010297</a> , <a href="#">BC072594</a> , <a href="#">BC172737</a> , <a href="#">BU705617</a>
Synonyms:	AW558298; CBP; CBP/p300; KAT3A; p300/CBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)