

Product datasheet for **MP203003**

Cryge Mouse qPCR Primer Pair (NM_007777)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12968
Forward Sequence:	ACCCTGACTACCAGCAGTGGAT
Reverse Sequence:	AGCAGTCGTCTGTGATCTCCAC
ACCN:	BC078417 , BC132506 , BC132508 , NM_007777 , NM_007777.1 , NM_007777.2 , NM_007777.3 , BC057012
UniProt ID:	Q03740
Synonyms:	Cryg-2; Cryg-6; DGcry-2; Elo
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

