

Product datasheet for **MP203023**

Cttnbp2nl Mouse qPCR Primer Pair (NM_030249)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80281
Forward Sequence:	CTGTCCTCAAGCGGTGATGAA
Reverse Sequence:	TCTCCTTCAGCAGTGTCTGTG
ACCN:	<u>BC003236</u> , <u>NM_030249</u> , <u>NM_030249.1</u> , <u>NM_030249.2</u> , <u>NM_030249.3</u> , <u>NM_030249.4</u> , <u>BC006952</u> , <u>BY195820</u>
UniProt ID:	<u>Q99LJO</u>
Synonyms:	AA552995; AA589392; AU018624; BC003236; mKIAA1433
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

