

## Product datasheet for **MP203038**

### Clec4a2 Mouse qPCR Primer Pair (NM\_011999)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 26888                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | CAAGAGTGAGGAGAACTGCTCC                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | GCAGCATGAGTGTCCAAGATCC                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <u><a href="#">BC006623</a></u> , <u><a href="#">BC075729</a></u> , <u><a href="#">BC096565</a></u> , <u><a href="#">NM_011999</a></u> , <u><a href="#">NM_011999.1</a></u> , <u><a href="#">NM_011999.2</a></u> , <u><a href="#">NM_011999.3</a></u> ,<br><u><a href="#">NM_011999.4</a></u>                                                 |
| UniProt ID:       | <u><a href="#">Q9QZ15</a></u>                                                                                                                                                                                                                                                                                                                 |
| Synonyms:         | Clec4a; Clec4f6; DCIR; Dcir1                                                                                                                                                                                                                                                                                                                  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions).<br>Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |

