

Product datasheet for **MP203066**

Cul5 Mouse qPCR Primer Pair (NM_027807)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75717
Forward Sequence:	GTGTAAATGCGCTGGTGACCTC
Reverse Sequence:	CTCCAAGTCCTTCAACATCGGC
ACCN:	<u>BC075710</u> , <u>NM_027807</u> , <u>NM_027807.1</u> , <u>NM_027807.2</u> , <u>NM_027807.3</u> , <u>BC019661</u> , <u>BC030380</u> , <u>BC037054</u>
UniProt ID:	<u>Q9D5V5</u>
Synonyms:	4921514I20Rik; 8430423K24Rik; AI852817; C030032G03Rik; C330021I08Rik; VACM-1; VACM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)