

Product datasheet for **MP203067**

Cul7 Mouse qPCR Primer Pair (NM_025611)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66515
Forward Sequence:	GCAGGAACAGAATCCATGACTGG
Reverse Sequence:	GTCTTTCCAGGACGCAGTCAGA
ACCN:	<u>NM_025611</u> , <u>NM_025611.1</u> , <u>NM_025611.2</u> , <u>NM_025611.4</u> , <u>NM_025611.5</u> , <u>BC138998</u> , <u>BC019645</u> , <u>BC024345</u> , <u>BC026946</u> , <u>BC059865</u> , <u>BC065096</u> , <u>BC100544</u> , <u>BC139000</u> , <u>BY749504</u>
UniProt ID:	<u>Q8VE73</u>
Synonyms:	2510004L20Rik; AA409809; C230011P08Rik; p185; p193
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)