

Product datasheet for **MP203069**

Cuta Mouse qPCR Primer Pair (NM_026307)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67675
Forward Sequence:	GCTGTCAATGGCTTCTGGAAGC
Reverse Sequence:	TGCCAGACGCTTCTTACCACT
ACCN:	BC024422 , NM_026307 , NM_026307.1 , NM_026307.2 , NM_026307.3 , BC021463 , BC028531 , BC037500 , BY033979 , BY057392
UniProt ID:	Q9CQ89
Synonyms:	0610039D01Rik; 1810022E02Rik; 1810060C03Rik; 2700094G22Rik; AI326454
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)