

Product datasheet for **MP203071**

Cux1 Mouse qPCR Primer Pair (NM_009986)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13047
Forward Sequence:	CCAAGAACGGCATTGCCCAGAG
Reverse Sequence:	CCTTTCTGGGTCAGTTTGCTCC
ACCN:	BC090847 , NM_009986 , NM_009986.1 , NM_009986.2 , NM_009986.3 , NM_009986.4 , BC014289 , BC054370 , BC079570 , BM944198 , BP762606 , BX522962
UniProt ID:	P53564
Synonyms:	CDP; Cutl1; Cux; Cux-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)