

## Product datasheet for **MP203091**

### Dmtf1 Mouse qPCR Primer Pair (NM\_011806)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23857
Forward Sequence:	ATGACGGCAACTACAGAGGTGG
Reverse Sequence:	ACGCTTGGCTAACTGCTGAGAC
ACCN:	<u><a href="#">BC045141</a></u> , <u><a href="#">NM_011806</a></u> , <u><a href="#">NM_011806.1</a></u> , <u><a href="#">NM_011806.2</a></u> , <u><a href="#">NM_011806.3</a></u> , <u><a href="#">BC010272</a></u> , <u><a href="#">BC060109</a></u> , <u><a href="#">BF466048</a></u> , <u><a href="#">BM243709</a></u> , <u><a href="#">BY247901</a></u>
UniProt ID:	<u><a href="#">Q8CE22</a></u>
Synonyms:	Dimp; Dmp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)