

Product datasheet for **MP203107**

Ccn1l Mouse qPCR Primer Pair (NM_019937)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56706
Forward Sequence:	CGTTCCAGATCACACAGCCCAA
Reverse Sequence:	GGCTTTCACGTGTGACTGCGAGA
ACCN:	<u>BC094383</u> , <u>NM_019937</u> , <u>NM_019937.1</u> , <u>NM_019937.2</u> , <u>NM_019937.3</u> , <u>BC007177</u> , <u>BC037204</u> , <u>BC069909</u> , <u>BC117816</u> , <u>NM_019937.4</u>
UniProt ID:	<u>Q52KE7</u>
Synonyms:	2610030E23Rik; ania-6a; AU018493; Ccn1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

