

Product datasheet for **MP203110**

Cnnm2 Mouse qPCR Primer Pair (NM_033569)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94219
Forward Sequence:	TAAGCAGACGGACAGCGAGATG
Reverse Sequence:	GACATCTGGGATGGGCTAAACG
ACCN:	<u>BC125033</u> , <u>NM_033569</u> , <u>NM_033569.1</u> , <u>NM_033569.2</u> , <u>NM_033569.3</u> , <u>BC012956</u> , <u>BC027387</u> , <u>BC038027</u> , <u>BC052513</u>
UniProt ID:	<u>Q3TWN3</u>
Synonyms:	Acdp2; AU015877; AW048635; Clp2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

