

Product datasheet for **MP203160**

Cyhr1 Mouse qPCR Primer Pair (NM_019396)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54151
Forward Sequence:	GCACCAGAAAGAGGAATGCCAG
Reverse Sequence:	TCGTCGAGGATCTCCATCAGCT
ACCN:	<u>BC016196</u> , <u>BC025112</u> , <u>BC062184</u> , <u>BC082772</u> , <u>NM_019396</u> , <u>NM_019396.1</u> , <u>NM_019396.2</u> , <u>NM_019396.3</u> , <u>BC062184.1</u> , <u>BB849289</u> , <u>BC039646</u> , <u>BE653118</u> , <u>BE859486</u> , <u>BG087565</u> , <u>BM939179</u> , <u>BX512375</u> , <u>BX522828</u> , <u>BY098544</u> , <u>BY287678</u>
UniProt ID:	<u>Q9QXA1</u>
Synonyms:	1110031M01Rik; AU042374; Chrp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

