

Product datasheet for **MP203166**

Kcp Mouse qPCR Primer Pair (NM_001029985)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	333088
Forward Sequence:	TGCTTGACCAGGTCACAGAAC
Reverse Sequence:	GTGACACATGCAGGAGGAACAG
ACCN:	<u>NM_001029985</u> , <u>NM_001029985.1</u> , <u>NM_001029985.2</u> , <u>NM_001029985.3</u> , <u>NM_001029985.4</u> , <u>BC032979</u> , <u>BC059023</u> , <u>BC059909</u> , <u>BC119794</u> , <u>BQ179923</u>
UniProt ID:	<u>Q3U492</u>
Synonyms:	AW060220; CRIM-2; Crim2; Gm793; KCP-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

