

Product datasheet for **MP203179**

Cypt1 Mouse qPCR Primer Pair (NM_025738)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66742
Forward Sequence:	CGTCGTTTTCTAGCAGAACTAC
Reverse Sequence:	GTTCGATCAGCTCCAGTTGGCT
ACCN:	BC061079 , NM_025738 , NM_025738.1 , NM_025738.2 , NM_025738.3 , BC048455 , BU961231
UniProt ID:	Q8CH20
Synonyms:	4931420D14Rik; Ckt1r3; Cypt5; Cypt6; TISP55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

