

Product datasheet for **MP203187**

Crisp3 Mouse qPCR Primer Pair (NM_009639)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11572
Forward Sequence:	TGTGGAGTTGCTGAATGCCCTG
Reverse Sequence:	ATGGTGCTCTTGCTGTGTAAGCT
ACCN:	<u>BC022573</u> , <u>BC132536</u> , <u>BC132538</u> , <u>NM_009639</u> , <u>NM_009639.1</u> , <u>NM_009639.2</u> , <u>BC114345</u> , <u>NM_009639.3</u>
UniProt ID:	<u>Q03402</u>
Synonyms:	Aeg2; CRISP-3; CRS3; SGP28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

