

Product datasheet for **MP203196**

Cysltr2 Mouse qPCR Primer Pair (NM_133720)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70086
Forward Sequence:	CATGCTTGGAGCTGAGTCCTCA
Reverse Sequence:	TAGCAGACGGTGAGCGTGAGAA
ACCN:	<u>BC132414, BC132418, NM_133720, NM_133720.1, NM_133720.2, NM_133720.3, BE848050</u>
UniProt ID:	<u>Q920A1</u>
Synonyms:	2300001H05Rik; Cltr2; Cyslt2; CYSLT2R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

