

Product datasheet for **MP203211**

Cyba Mouse qPCR Primer Pair (NM_007806)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13057
Forward Sequence:	GTCATCTGTCTGCTGGAGTATC
Reverse Sequence:	CGGACGTAGTAATTCCTGGTGAG
ACCN:	<u>BC026791</u> , <u>NM_007806</u> , <u>NM_007806.1</u> , <u>NM_007806.2</u> , <u>NM_007806.3</u> , <u>BC060635</u> , <u>BY005999</u>
UniProt ID:	<u>Q61462</u>
Synonyms:	b558; nmf333; p22-phox; p22phox
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)