

Product datasheet for **MP203218**

Cyb5r3 Mouse qPCR Primer Pair (NM_029787)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109754
Forward Sequence:	CGTTGTCAGGACGGTGAAGTCT
Reverse Sequence:	TAGCACACGGTGTGGTCGTTTG
ACCN:	<u>BC004760</u> , <u>BC032013</u> , <u>BC043074</u> , <u>NM_029787</u> , <u>NM_029787.1</u> , <u>NM_029787.2</u> , <u>BC018574</u> , <u>NM_029787.3</u>
UniProt ID:	<u>Q9DCN2</u>
Synonyms:	0610016L08Rik; 2500002N19Rik; C85115; Dia-1; Dia1; WU:AL591952.1-001; WU:AL591952.1-002
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)