

## Product datasheet for **MP203220**

### Cyb5b Mouse qPCR Primer Pair (NM\_025558)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 66427                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | AAGCGAAACTCGGCGGAAGAGA                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | TGTTCCAGCAGAACCTCTTCCC                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <u>BC054749, BC058812, BC062980, NM_025558, NM_025558.1, NM_025558.2, NM_025558.3, NM_025558.4, NM_025558.5, BC046479, BC046810, BC048147</u>                                                                                                                                                                                                 |
| UniProt ID:       | <u>Q9CQX2</u>                                                                                                                                                                                                                                                                                                                                 |
| Synonyms:         | 1810044O22Rik; AU015618; Cyb5m                                                                                                                                                                                                                                                                                                                |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |

