

Product datasheet for **MP203263**

Cyp2c29 Mouse qPCR Primer Pair (NM_007815)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13095
Forward Sequence:	GCTCTCCTACTCCTGCTGAAGT
Reverse Sequence:	ATGTGGCTCCTGTCTTGCATGC
ACCN:	<u>BC013895</u> , <u>BC019908</u> , <u>NM_007815</u> , <u>NM_007815.1</u> , <u>NM_007815.2</u> , <u>NM_007815.3</u>
UniProt ID:	<u>Q64458</u>
Synonyms:	Ah-2; Ahh-1; AHOH; AHOHase; AII59681; AII96010; Cyp2c; P450-2C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

