

Product datasheet for **MP203267**

Cyp2c23 Mouse qPCR Primer Pair (NM_001001446)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226143
Forward Sequence:	GGCTATGATGTGGTGAAGGAAGC
Reverse Sequence:	TCGCAGGAGTTTCCACCTCTCT
ACCN:	BC025819 , BC026766 , BC034834 , NM_001001446 , NM_001001446.1 , NM_001001446.2 , NM_001001446.3 , BC022141 , BC024068 , BK005218
UniProt ID:	A0A1D5RLI1
Synonyms:	AW107714; b2b386Clo; BC034834; Cyp2c44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

