

Product datasheet for **MP203279**

Cyp2d22 Mouse qPCR Primer Pair (NM_019823)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56448
Forward Sequence:	CCGAGATTTGACTGATGCCTTCC
Reverse Sequence:	GCAGAGAACAGGTCACCCACTA
ACCN:	<u>BC016256</u> , <u>NM_019823</u> , <u>NM_019823.1</u> , <u>NM_019823.2</u> , <u>NM_019823.3</u> , <u>NM_019823.4</u> , <u>BC006853</u> , <u>BC025224</u> , <u>BX520002</u> , <u>BY297294</u>
Synonyms:	2D22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

