

Product datasheet for **MP203325**

Cyp5l Mouse qPCR Primer Pair (NM_020010)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13121
Forward Sequence:	ATCCAGAAGCGCAGGCTGTCAA
Reverse Sequence:	CAGTCCGATGAGCATCCCTGAT
ACCN:	BC031813 , NM_020010 , NM_020010.1 , NM_020010.2
UniProt ID:	Q8K0C4
Synonyms:	AI426508; Cyp51a1; CYPLI; Ldm; P450-14DM; P450LI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

