

Product datasheet for **MP203333**

Crlf1 Mouse qPCR Primer Pair (NM_018827)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12931
Forward Sequence:	CCTCTTCCAAGCCAAGTACCAG
Reverse Sequence:	CTTGGACGAAGTAAACGGTGCC
ACCN:	<u>NM_018827</u> , <u>NM_018827.2</u> , <u>NM_018827.3</u> , <u>BC107376</u> , <u>BC107377</u> , <u>BC114407</u> , <u>BC160203</u> , <u>NM_018827.4</u>
UniProt ID:	<u>Q9JIM58</u>
Synonyms:	CLF-; CLF-1; CRL; CRLM-3; CRLM3; NR6; NR6.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)