

Product datasheet for **MP203348**

Ckap4 Mouse qPCR Primer Pair (NM_175451)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216197
Forward Sequence:	GCTGGTATGTCCATCACGTCCT
Reverse Sequence:	CGAAAGTCCCAAATGTGGCTTGC
ACCN:	<u>NM_175451</u> , <u>NM_175451.1</u> , <u>BB033173</u> , <u>BB620334</u> , <u>BC025522</u> , <u>BC138318</u> , <u>BC138320</u> , <u>BE647088</u>
UniProt ID:	<u>Q8BMK4</u>
Synonyms:	5630400A09Rik; CLIMP-63; P63
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)