

Product datasheet for **MP203358**

Dpf3 Mouse qPCR Primer Pair (NM_058212)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70127
Forward Sequence:	CCTCTCAGGAAGACCACGACAA
Reverse Sequence:	CCAGGTGAGTATGACGCTAGTG
ACCN:	BC048572 , NM_058212 , NM_058212.1 , NM_058212.2 , BC048572.1 , BF469100
UniProt ID:	P58269
Synonyms:	2810403B03Rik; 6530402L11Rik; BAF45C; C78788; cer-d4; CERD4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

