

Product datasheet for **MP203363**

Ddb2 Mouse qPCR Primer Pair (NM_028119)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107986
Forward Sequence:	GCAAAAGTCGCAGAGGTCCTCA
Reverse Sequence:	CCTGACAATGCTGCTACTCCGT
ACCN:	BC099410 , NM_028119 , NM_028119.2 , NM_028119.3 , NM_028119.4 , NM_028119.5 , BC028955 , BC058589 , BC080821 , BU515611
UniProt ID:	Q99J79
Synonyms:	2610043A19Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)